

Clinical Pathology - Core Laboratory Room 3702, University Hospital 750 E. Adams Street Syracuse, NY 13210

#### **MEMORANDUM**

Kul Bol

TO: All Physicians

FROM: Katalin Banki, M.D., Director of Core Laboratory

Stephen Gwilt, MS,MT(ASCP)SH,SC, Supervisor Core Laboratory

DATE: July 24, 2017

RE: Changes to Hemoglobin A1c, C-reactive Protein, and Lactate

\_\_\_\_\_

### 1. New Method to Measure Hemoglobin A1c:

Old Method: HPLC

New Method: Immunoassay (Turbidimetric Inhibition Immunoassay) EPIC Code:

**LAB90** 

#### What is staying the same:

- Test results, reference values, and therapeutic goals, as both methods are traceable to DCCT/NGSP standards.
- HbA1c is lower in hemolytic anemias and increased in polycythemia and postsplenectomy.
- Can be used in Hemoglobin S, C, E, D traits.

#### What is changing:

- Performed: Everyday
- Acceptable Samples: EDTA-Lavender top tubes
- >10% HbF interfere with the assay

### 2. <u>Two Tests to Measure C-reactive Protein</u>:

- C-reactive protein in inflammatory conditions (CRP), Test Code CRP1
- Cardiac C-reactive protein, high-sensitivity in cardiac risk assessment (hs-CRP) LAB150
  - Serum/plasma samples are tested at different dilutions in the two tests, in order to reach an ideal precision and fast turn-around-time.
  - The test principle is the same (particle enhanced immunoturbidimetric method).

# Page 2

## Changes to Hemoglobin A1c, C-reactive Protein, and Lactate

- $\circ$  CRP <5 mg/L
- o hs-CRP <3 mg/L (<1.0 Low Risk; 1.0-3.0 Average Risk; >3.0 High Risk)

# 3. <u>Lactate: New Reagent with Increased Sensitivity:</u>

- Lactate levels will be estimated 15-20% higher with our new reagent, when lactate is >5 mmol/L.
- There will be only a modest increase (<5%) in the measured value, when lactate is <5 mmol/L.
- Consequently, the normal range of 0.5-2.2 mmol/L, and the critical value of 3.9 mmol/L are not changing.