

---

<b>Issue Date:</b>	<b>November 17, 1998</b>	<b>Procedure Number:</b> Evs136
<b>Last Revision Date:</b>	<b>July 14, 2001</b>	<b>Approved by:</b> Department Director
<b>Last Review Date:</b>	<b>March 27, 2008</b>	<b>Page:</b> 1 of 1

---

## **Restrooms: Toilet Bowl Cleaning**

**Purpose:** To keep toilet bowls within the facility, free of dirt, stains and body fluids by cleaning them often with disinfectant cleaner.

**Procedure:**

1. Mix quat-stat solution in accordance with manufacturer's instructions.
2. Pour the quat-stat solution from the one gallon pail into the toilet.
3. Using the bowl swab, push water from bowl into the drain. The lower water levels aid in cleaning the entire inside of the bowl.
4. With the swab, clean inside the bowl, starting at the lip. Look for mineral deposits and other stains. If they cannot be removed by the bowl swab, using quat-stat solution, ask your supervisor for toilet bowl cleaner.
5. Prepare a clean pail of quat-stat. Use clean rag, wrung out in quat-stat solution, to clean the remainder of the bowl; toilet seat and lid, top, bottom, and outside of bowl.
6. Once the toilet is cleaned, the remaining quat-stat cleaner should be placed in the toilet and flushed. The rag(s) should be placed in a soiled linen bag and placed in the appropriate container for soiled rags and mops.
7. **Sign off on cleaning sheet if applicable.**

**References:** None Cited